


SIALKOT · PAKISTAN

Manufacturer & Exporter

MABAN SURGICAL

Surgical Instruments Catalogue

Scissors · Forceps · Clamps · Retractors · Needle Holders




info@mabansurgical.online · WhatsApp +92 345 6776994

mabansurgical.online

Surgical Instruments

Scissors · Forceps · Clamps · Retractors · Needle Holders



In This Section

- | | |
|--|--|
| 01 Mayo Dissecting Scissors — Straight | 14 Crile Hemostatic Forceps |
| 02 Metzenbaum Dissecting Scissors — Curved | 15 Babcock Intestinal Forceps |
| 03 Iris Scissors — Straight / Curved | 16 DeBakey Atraumatic Tissue Forceps |
| 04 Halsted Mosquito Forceps — Straight / Curved | 17 Russian Tissue Forceps |
| 05 Kelly Hemostatic Forceps — Straight / Curved | 18 Olsen-Hegar Needle Holder |
| 06 Allis Tissue Forceps | 19 Derf Needle Holder |
| 07 Backhaus Towel Clamp | 20 Senn-Miller Retractor |
| 08 Mayo-Hegar Needle Holder | 21 Volkman Rake Retractor |
| 09 Crile-Wood Needle Holder | 22 Weitlaner Self-Retaining Retractor |
| 10 Army-Navy Retractor Set | 23 Scalpel Handle — No. 3 / No. 4 |
| 11 Adson Tissue Forceps — Toothed | 24 Yankauer Suction Tube |
| 12 Spencer Stitch Scissors | 25 Littauer Stitch Scissors |
| 13 Rochester-Pean Hemostatic Forceps | |

Mayo Dissecting Scissors — Straight

Code: MS-SUR-001 | Scissors | General Surgery

Heavy-pattern dissecting scissors for cutting dense tissue and suture material. Straight pattern with blunt-blunt blades.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	14 cm, 17 cm, 19 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Metzenbaum Dissecting Scissors — Curved

Code: MS-SUR-002 | Scissors | General Surgery

Delicate-pattern scissors for cutting fine tissue. Curved blades with blunt tips for controlled dissection.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	12 cm, 14 cm, 18 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Iris Scissors — Straight / Curved

Code: MS-SUR-003 | Scissors | General & Specialty Surgery

Fine, precision scissors originally for ophthalmic use, widely used for delicate dissection in multiple specialties.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	9 cm, 11.5 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Halsted Mosquito Forceps — Straight / Curved

Code: MS-SUR-004 | Forceps | General Surgery

Small hemostatic forceps for clamping fine vessels. Fully serrated jaws with precision box lock.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	12.5 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Kelly Hemostatic Forceps — Straight / Curved

Code: MS-SUR-005 | Forceps | General Surgery

Medium-size hemostats for clamping larger vessels and tissue. Transverse jaw serrations.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	14 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Allis Tissue Forceps

Code: MS-SUR-006 | Forceps | General Surgery

Grasping forceps with interlocking teeth for holding dense tissue such as fascia and breast tissue.




Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	15 cm, 19 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Backhaus Towel Clamp

Code: **MS-SUR-007** | **Clamps** | **General Surgery**

Sharp-pointed towel clamp for securing drapes and towels. Curved jaws with positive locking ratchet.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	11 cm, 13 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Mayo-Hegar Needle Holder

Code: **MS-SUR-008** | **Needle Holders** | **General Surgery**

Heavy-pattern needle holder for closing dense tissue. Cross-serrated jaws with hinged box lock.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	15 cm, 18 cm, 20 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Crile-Wood Needle Holder

Code: **MS-SUR-009** | **Needle Holders** | **General Surgery**

Delicate needle holder for fine suture work. Slim jaws with precision serrations.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	15 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Army-Navy Retractor Set

Code: **MS-SUR-010** | **Retractors** | **General Surgery**

Double-ended handheld retractors for superficial wound retraction. Supplied as a pair.



Material	Stainless Steel
Finish	Satin Finish
Sizes	21 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Adson Tissue Forceps — Toothed

Code: MS-SUR-011 | Forceps | General Surgery

Thumb forceps with fine 1×2 teeth for grasping delicate tissue during closure.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	12 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Spencer Stitch Scissors

Code: MS-SUR-012 | Scissors | Suture Removal

Short scissors with a notched lower blade designed to slide under sutures for safe removal.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	9 cm, 11 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Rochester-Pean Hemostatic Forceps

Code: MS-SUR-013 | Forceps | General Surgery

Large hemostatic forceps with full-length jaw serrations for clamping heavy tissue and vessels.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	16 cm, 18 cm, 20 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Crile Hemostatic Forceps

Code: MS-SUR-014 | Forceps | General Surgery

Medium hemostats with fully serrated jaws for general clamping of vessels and tissue.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	14 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Babcock Intestinal Forceps

Code: MS-SUR-015 | Forceps | General & Abdominal Surgery

Fenestrated atraumatic forceps for grasping tubular structures such as bowel and fallopian tubes.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	16 cm, 20 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

DeBakey Atraumatic Tissue Forceps

Code: MS-SUR-016 | Forceps | Vascular & General Surgery

Thumb forceps with fine atraumatic serrations for delicate vascular and soft-tissue handling.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	15 cm, 19 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Russian Tissue Forceps

Code: **MS-SUR-017** | **Forceps** | **General Surgery**

Thumb forceps with a broad serrated tip for a secure grip on dense tissue and dressings.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	15 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Olsen-Hegar Needle Holder

Code: **MS-SUR-018** | **Needle Holders** | **General Surgery**

Combination needle holder and scissors for efficient suturing without changing instruments.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	14 cm, 17 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Derf Needle Holder

Code: MS-SUR-019 | Needle Holders | General & Plastic Surgery

Short, delicate needle holder for fine suture work in confined operative areas.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	12 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Senn-Miller Retractor

Code: MS-SUR-020 | Retractors | General & Small-Incision Surgery

Double-ended retractor with an L-shaped blade and rake prongs for small wound retraction.



Material	Stainless Steel
Finish	Satin Finish
Sizes	16 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Volkman Rake Retractor

Code: **MS-SUR-021** | Retractors | Orthopedic & General Surgery

Pronged rake retractor for firm tissue retraction. Available with blunt or sharp prongs.



Material	Stainless Steel
Finish	Satin Finish
Sizes	21 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Weitlaner Self-Retaining Retractor

Code: **MS-SUR-022** | Retractors | General & Orthopedic Surgery

Self-retaining wound retractor with hinged arms, prongs and an adjustable ratchet lock.



Material	Stainless Steel
Finish	Satin Finish
Sizes	13 cm, 16 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Scalpel Handle — No. 3 / No. 4

Code: MS-SUR-023 | Accessories | General Surgery

Reusable scalpel handle accepting standard detachable surgical blades. Graduated scale on flat pattern.



Material	Stainless Steel
Finish	Satin Finish
Sizes	No. 3, No. 4
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Yankauer Suction Tube

Code: MS-SUR-024 | Accessories | Suction & Aspiration

Rigid suction tube with a large-bore tip for oropharyngeal suctioning and fluid removal.



Material	Stainless Steel
Finish	Satin Finish
Sizes	Standard
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

Littauer Stitch Scissors

Code: MS-SUR-025 | Scissors | Suture Removal

Lightweight stitch scissors with a long thin blade for reaching and removing sutures.



Material	Stainless Steel
Finish	Mirror Polish, Satin Finish
Sizes	14 cm
Custom	Sizes, finishes, marking, packaging
Packaging	Standard / Individual / Custom / Private-label
MOQ	Confirm in RFQ
Lead time	Confirm in RFQ

◆ CUSTOMIZABLE ◆ OEM / ODM

START SOURCING

Have a Product Requirement?

Tell us what you need. We can discuss specifications, customization, branding, packaging and bulk supply requirements.

ADDRESS

Laalpur Head, Marala Road, Sialkot, Pakistan

EMAIL

info@mabansurgical.online

WHATSAPP

+92 345 6776994

WEB

mabansurgical.online



Visit Website



WhatsApp Us